


DOUBLE MOUNTING PIN



| CODE | A | B | C | D | E | F | G | H | I | J |
|-------------|------|------|----|-----|-----|----|----|-----|-----|-----|
| MEDMP 07003 | 25 | 19 | 40 | 70 | 100 | 80 | 20 | 12 | 12 | 210 |
| MEDMP 07005 | 25 | 19 | 38 | 96 | 78 | 62 | 22 | 12 | 12 | 212 |
| MEDMP 07006 | 25 | 19 | 26 | 66 | 78 | 70 | 18 | 6.5 | 6.5 | 170 |
| MEDMP 07008 | 28 | 22 | 35 | 75 | 90 | 75 | 20 | 12 | 12 | 200 |
| MEDMP 07009 | 28 | 22 | 40 | 70 | 100 | 80 | 20 | 12 | 12 | 210 |
| MEDMP 07010 | 28 | 22 | 40 | 90 | 100 | 80 | 20 | 12 | 12 | 230 |
| MEDMP 07012 | 28.6 | 22.2 | 40 | 89 | 100 | 80 | 20 | 12 | 12 | 22 |
| MEDMP 07014 | 32 | 25 | 40 | 115 | 95 | 81 | 26 | 12 | 12 | 250 |
| MEDMP 07016 | 36.5 | 28 | 36 | 150 | 91 | 81 | 26 | 12 | 12 | 277 |
| MEDMP 07018 | 36.5 | 28 | 45 | 150 | 95 | 80 | 28 | 12 | 12 | 290 |